**3M**

**Cavilon™**

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**A prospective randomized trial of the effect of a soluble adhesive on the ease of dressing removal following hypospadias repair** pag. 8

Sanders C, Young A, McAndrew HF, Kenny SE. A prospective randomized trial of the effect of a soluble adhesive on the ease of dressing removal following hypospadias repair. *J Pediatr Urol.* 2007 Jun;3(3):209-213. doi:10.1016/j.jpuro.2006.08.006. Epub 2006 Dec 12. PMID: 18947737.

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**Applying skin barrier film for skin tear management in patients with central venous catheterization** **V<sub>A</sub>** pag. 9

Chen YH, Hsieh HL, Shih WM. Applying skin barrier film for skin tear management in patients with central venous catheterization. *Adv Skin Wound Care.* 2020 Nov;33(11):582-586. doi:10.1097/01.ASW.0000717208.20481.a0. PMID: 33065679.

**Effect of an acrylic terpolymer barrier film beneath transparent catheter dressings on skin integrity, risk of dressing disruption, catheter colonisation and infection** **V<sub>A</sub>** pag. 10

Pivkina AI, Gusarov VG, Blot SI, Zhivotneva IV, Pasko NV, Zamyatin MN. Effect of an acrylic terpolymer barrier film beneath transparent catheter dressings on skin integrity, risk of dressing disruption, catheter colonisation and infection. *Intensive Crit Care Nurs.* 2018 Jun;46:17-23. doi:10.1016/j.iccn.2017.11.002. Epub 2018 Mar 23. PMID: 29576395.

**Use of a barrier film (3M Cavilon No Sting Barrier Film) to reduce local skin complications around peripherally inserted central catheter lines: a randomised prospective controlled study** **V<sub>A</sub>** pag. 11

George M, Pal U, Guduri V, et al. Use of a barrier film (3M Cavilon No Sting Barrier Film) to reduce local skin complications around peripherally inserted central catheter lines: a randomised prospective controlled study. *WCET Journal.* 2016;36(4):8-13.

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**Skin impairment associated with vascular access devices and semi-permeable transparent dressings** **V<sub>A</sub>** pag. 12

Hitchcock J, Anderson L, Escorcio J, et al. Skin impairment associated with vascular access devices and semi-permeable transparent dressings. Imperial College Healthcare NHS Trust. *NIVAS*. 2015. Poster Presentation.

**In search of a better central line dressing protocol in the autologous bone marrow reinfusion patient** **V<sub>A</sub>** pag. 13

Link D, Cutler C. In search of a better central line dressing protocol in the autologous bone marrow reinfusion patient. 3M Clinical study (1998). White Paper.

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**An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence: economics of skin damage prevention** **I** pag. 16

Bliss DZ, Zehrer C, Savik K, Smith G, Hedblom E. An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence: economics of skin damage prevention. *J Wound Ostomy Continence Nurs*. 2007 Mar-Apr;34(2):143-152; discussion 152. doi:10.1097/01.WON.0000264825.03485.40. PMID: 17413828.

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### **PERILESIONALE**

**Effectiveness of association of multilayer compression therapy and periwound protection with Cavilon™ (No Sting Barrier Film) in the treatment of venous leg ulcers** **P<sub>W</sub>** pag. 18

Serra N, Palomar-Llatas F, Pujalte BF, et al. Effectiveness of association of multilayer compression therapy and periwound protection with Cavilon™ (No Sting Barrier Film) in the treatment of venous leg ulcers. *Gerokomos*. 2010;21(3):124-130. [http://scielo.isciii.es/scielo.php?pid=S1134-928X2010000300006&script=sci\\_abstract&tlng=en](http://scielo.isciii.es/scielo.php?pid=S1134-928X2010000300006&script=sci_abstract&tlng=en).

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### The protective effects of a new preparation on wound edges

pag. 19

Neander KD, Hesse F. The protective effects of a new preparation on wound edges. *J Wound Care*. 2003 Nov;12(10):369-371. doi:10.12968/jowc.2003.12.10.26548

### Comparison of two peri-wound skin protectants in venous leg ulcers: a randomised controlled trial

pag. 20

Cameron J, Hoffman D, Wilson J, Cherry G. Comparison of two peri-wound skin protectants in venous leg ulcers: a randomised controlled trial. *J Wound Care*. 2005 May;14(5):233-236. doi:10.12968/jowc.2005.14.5.26779

### A liquid film-forming acrylate for peri-wound protection: a systematic review and meta-analysis (3M Cavilon No Sting Barrier Film)

pag. 21

Schuren J, Becker A, Sibbald RG. A liquid film-forming acrylate for peri-wound protection: a systematic review and meta-analysis (3M Cavilon No Sting Barrier Film). *Int Wound J*. 2005 Sep;2(3):230-238. doi:10.1111/j.1742-4801.2005.00131.x

### Comparative study of a barrier product versus zinc oxide for the treatment of incontinent lesions

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### Comparing cost per use of 3M Cavilon No Sting Barrier Film with zinc oxide oil in incontinent patients

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#### Randomized, paired comparison of No Sting Barrier Film versus sorbolene cream (10% glycerine) skin care during postmastectomy irradiation

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### **Randomized control trial of 3M™ Cavilon™ No Sting Barrier Film for the prevention of radiation dermatitis in patients with nasopharyngeal carcinoma**

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## **PRESTAZIONI**

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**pag. 28**

### **A comparison of the durability of four barrier film products over a 72 hour period on human volunteers**

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### **A clinical evaluation of 3M No Sting Barrier Film**

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2. Dati 3M in archivio. CLIN-RPT-FINAL-ICH2-US-05-291160.

3. Dati 3M in archivio. TEAM-MISC-05-005732 and SPONSOR FINAL RPT-05-002049.

4. Campbell K, Woodbury MG, Whittle H, Labate T, Hoskin A. A clinical evaluation of 3M Cavilon No Sting Barrier Film. Ostomy Wound Manage. 2000;46(1)24-30.

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

5. McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. *J Wound Ostomy Continence Nurs.* 2013;40(4):365-380.

6. Dati 3M in archivio. CLIN-MISC-US-05-202211.

# A prospective randomized trial of the effect of a soluble adhesive on the ease of dressing removal following hypospadias repair

Sanders C, Young A, McAndrew HF, Kenny SE. A prospective randomized trial of the effect of a soluble adhesive on the ease of dressing removal following hypospadias repair. *J Pediatr Urol.* 2007 Jun;3(3):209-213. doi:10.1016/j.jpuro.2006.08.006. Epub 2006 Dec 12. PMID: 18947737.

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
# Applying skin barrier film for skin tear management in patients with central venous catheterization.

Chen YH, Hsieh HL, Shih WM. Applying skin barrier film for skin tear management in patients with central venous catheterization. *Adv Skin Wound Care*. 2020 Nov;33(11):582-586. doi:10.1097/01.ASW.0000717208.20481.a0. PMID: 33065679.

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
# Effect of an acrylic terpolymer barrier film beneath transparent catheter dressings on skin integrity, risk of dressing disruption, catheter colonisation and infection

Pivkina AI, Gusarov VG, Blot SI, Zhivotneva IV, Pasko NV, Zamyatin MN. Effect of an acrylic terpolymer barrier film beneath transparent catheter dressings on skin integrity, risk of dressing disruption, catheter colonisation and infection. *Intensive Crit Care Nurs*. 2018 Jun;46:17-23. doi:10.1016/j.iccn.2017.11.002. Epub 2018 Mar 23. PMID: 29576395.

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# Use of a barrier film (3M Cavilon No Sting Barrier Film) to reduce local skin complications around peripherally inserted central catheter lines: a randomised prospective controlled study

George M, Pal U, Guduri V, et al. Use of a barrier film (3M Cavilon No Sting Barrier Film) to reduce local skin complications around peripherally inserted central catheter lines: a randomised prospective controlled study. *WCET Journal*. 2016;36(4):8-13.

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
## Incidenze minori

# Skin impairment associated with vascular access devices and semi-permeable transparent dressings

Hitchcock J, Anderson L, Escorcio J, et al. Skin impairment associated with vascular access devices and semi-permeable transparent dressings. Imperial College Healthcare NHS Trust. NIVAS. 2015. Poster Presentation.

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# In search of a better central line dressing protocol in the autologous bone marrow reinfusion patient

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# An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence: economics of skin damage prevention

Bliss DZ, Zehrer C, Savik K, Smith G, Hedblom E. An economic evaluation of four skin damage prevention regimens in nursing home residents with incontinence: economics of skin damage prevention. *J Wound Ostomy Continence Nurs.* 2007 Mar-Apr;34(2):143-152; discussion 152. doi:10.1097/01.WON.0000264825.03485.40. PMID: 17413828.

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
AUTORE	TITOLO
<b>Beeckman</b>	Prevention and treatment of incontinence-associated dermatitis: literature review
<b>Bernatchez SF</b>	Reducing friction on skin at risk: the use of 3M™ Cavilon™ No Sting Barrier Film
<b>Bonnetblanc</b>	Periwound protection with 3M™ Cavilon™ No Sting Barrier Film for patients with chronic venous leg ulcers, a randomised multi-centre trial
<b>Coutts</b>	Periwound skin protection: a comparison of a new skin barrier vs. traditional therapies in wound management
<b>Gerber</b>	Wound protection for delicate baby and children's skin
<b>Guest JF</b>	Clinical and economic evidence supporting a transparent barrier film dressing in incontinence-associated dermatitis and periwound skin protection
<b>Guest JF</b>	Relative cost-effectiveness of a skin protectant in managing venous leg ulcers in the UK
<b>Hampton</b>	Film subjects win the day
<b>Hune</b>	Innovative strategy to prevent incontinence dermatitis
<b>Issberner</b>	A comparative study of the skin protectant performance of five barrier films
<b>LaVoie</b>	Comparison of the effectiveness of five different skin protective products
<b>Lopez JR</b>	Diaper rash. Local treatment with barrier products and quality of life
<b>Mooney</b>	Use of a non-alcohol incontinence barrier film on patients with severely compromised skin
<b>Murray</b>	Evaluation of routine use of an alcohol-free barrier film on patients with urinary and/or fecal incontinence
<b>Prather P</b>	Effectiveness of topical skin products in the treatment and prevention of incontinence-associated dermatitis: a systematic review
<b>Zehrer</b>	A comparison of cost and efficacy of three incontinence skin barrier products
<b>Zimmaro Bliss</b>	Incontinence-associated skin damage in nursing home residents: a secondary analysis of a prospective, multi-centre study

# Effectiveness of association of multilayer compression therapy and periwound protection with Cavilon™ (No Sting Barrier Film) in the treatment of venous leg ulcers

Serra N, Palomar-Llatas F, Pujalte BF, et al. Effectiveness of association of multilayer compression therapy and periwound protection with Cavilon™ (No Sting Barrier Film) in the treatment of venous leg ulcers. *Gerokomos*. 2010;21(3):124-130. [http://scielo.isciii.es/scielo.php?pid=S1134-928X2010000300006&script=sci\\_abstract&tlng=en](http://scielo.isciii.es/scielo.php?pid=S1134-928X2010000300006&script=sci_abstract&tlng=en).

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# A liquid film-forming acrylate for peri-wound protection: a systematic review and meta-analysis (3M Cavilon No Sting Barrier Film)

Schuren J, Becker A, Sibbald RG. A liquid film-forming acrylate for peri-wound protection: a systematic review and meta-analysis (3M Cavilon No Sting Barrier Film). *Int Wound J.* 2005 Sep;2(3):230-238. doi:10.1111/j.1742-4801.2005.00131.x

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# Comparative study of a barrier product versus zinc oxide for the treatment of incontinent lesions

Lopez JR, Perejamo MA, Torra JE, et al. Comparative study of a barrier product versus zinc oxide for the treatment of incontinent lesions. World Union of Wound Healing Societies (WUWHS) 2004. Oral Presentation.

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
# Comparing cost per use of 3M Cavilon No Sting Barrier Film with zinc oxide oil in incontinent patients

Baatenburg de Jong H, Admiraal H. Comparing cost per use of 3M Cavilon No Sting Barrier Film with zinc oxide oil in incontinent patients. *J Wound Care*. 2004 Oct;13(9):398-400. doi:10.12968/jowc.2004.13.9.27264. PMID: 15517755.

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


AUTORE	TITOLO
Bär	Ulcer edge protection with a polymer protective film
Bracelert	Clinical experience with an alcohol-free skin protectant 3M™ Cavilon™ No Sting Barrier Film
Campbell	The use of liquid film to treat severe incontinent dermatitis: case reports
Chan A	The use of a No Sting Barrier Film treatment protocol compared to routine clinical care for the treatment of stage 1 and 2 pressure injuries in long-term care
Garcia	3M Cavilon No Sting Barrier Film: an evaluation of periwounds prone to maceration
Gómez T	<i>In vivo</i> evaluation using confocal microscopy of protective effect of No Sting Barrier Film 3M Cavilon on periwound skin
Gonzalez	The use of 3M Cavilon No Sting Barrier Film to prevent maceration in pressure ulcers treated with an adhesive hydrocolloid dressing
Guest JF	Relative cost-effectiveness of a skin protectant in managing venous leg ulcers in the UK
Lazaro-Martinez JL	Reducing skin maceration in exudative diabetic foot ulcers

# Randomized, paired comparison of No Sting Barrier Film versus sorbolene cream (10% glycerine) skin care during postmastectomy irradiation

Graham P, Browne L, Capp A, et al. Randomized, paired comparison of No Sting Barrier Film versus sorbolene cream (10% glycerine) skin care during postmastectomy irradiation. *Int J Radiat Oncol Biol Phys.* 2004 Jan 1;58(1):241-246. doi:10.1016/s0360-3016(03)01431-7. PMID: 14697444.

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
# Randomized control trial of 3M™ Cavilon™ No Sting Barrier Film for the prevention of radiation dermatitis in patients with nasopharyngeal carcinoma

Chang L. Randomized control trial of 3M™ Cavilon™ No Sting Barrier Film for the prevention of radiation dermatitis in patients with nasopharyngeal carcinoma. Presented at: Third Congress of the World Union of Wound Healing Societies; June 2008; Toronto, Canada.

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
AUTORE	TITOLO
<b>Herst PM</b>	Protecting the radiation-damaged skin from friction: a mini review
<b>Kumar S</b>	Management of skin toxicity during radiation therapy: a review of the evidence
<b>Lam AC</b>	Phase III randomized pair comparison of a barrier film vs. standard skin care in preventing radiation dermatitis in post-lumpectomy patients with breast cancer receiving adjuvant radiation therapy
<b>Peskova</b>	Some uses of 3M Cavilon No Sting Barrier Film for the prevention of postradiation dermatitis in the head and neck region
<b>Shaw SZ</b>	Cavilon No Sting Barrier Film or topical corticosteroid (mometasone furoate) for protection against radiation dermatitis: a clinical trial

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# A comparison of the durability of four barrier film products over a 72 hour period on human volunteers

Houser T, Zerweck C, Grove G. A comparison of the durability of four barrier film products over a 72 hour period on human volunteers. Poster at Clinical Symposium for Advances in Skin and Wound Care (CSASWC), Orlando, Florida; 2010.

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
## A clinical evaluation of 3M No Sting Barrier Film

Campbell K, Woodbury MG, Whittle H, Labate T, Hoskin A. A clinical evaluation of 3M No Sting Barrier Film. *Ostomy Wound Manage.* 2000;46(1):24-30. PMID: 10732633.

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AUTORE	TITOLO
Bale S	The benefits of implementing a new skin care protocol in nursing homes
Garcia	Effectiveness of 3M Cavilon No Sting Barrier Film for preventing skin damage: a systematic review
Garrick V	A multidisciplinary team model of caring for patients with perianal Crohn's disease incorporating a literature review, topical therapy and personal practice
Grove	A motorized sliding sled apparatus for measuring the coefficient of friction of human skin <i>in vivo</i>
Hampton	The nursing care of common raw and bleeding conditions
Williams	3M Cavilon No Sting Barrier Film in the protection of vulnerable skin
Zehrer	Assessment of diaper-clogging potential of petrolatum moisture barriers

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© 2022 3M. Tutti i diritti riservati. 3M e